



### CENHBN-240ML, Floor Standing Type Milk Centrifuge

#### Features:

- Micro-computer control, LCD and digital dual display.
- Frequency conversion motor driving, maintenance-free
- Stainless steel centrifugal chamber, electric lid lock
- It can be operated by programming. Faults are diagnosed automatically.
- RCF can be directly set up and display, without RPM / RCF translated
- Imported well known frequency controller with precise speed control and durable in use
- Heating type optional.

**Code: ROTOR-A-8x30-1.4k**

#### Technical Parameters:

Model	CENHBN-240ML
Max Speed	1400 r/min
Max RCF	350 xg
Max Capacity	8X30ml angle rotor
Type Size	Ø25x215mm
Speed Accuracy	±10 r/min
Timer Range	0~99H59min/spin
Noise	≤60dB (A)
Power	AC 220V 50HZ 10A
Dimension	650 × 650 × 800 (L × W × H) mm
Weight	90 kg