



39240

39240, Waterproof Thermometer, -40°C to +200°C

Pocket-size with MAX/MIN recording.

- 2.75" (70mm) Stainless steel stem.
- °C/°F switchable ranges of -40 to 392°F (-40 to 200°C).
- MAX/MIN recall of high and low temperatures.
- Data Hold and Auto power Off features.
- Range: -40 to 392°F (-40 to 200°C).
Basic accuracy: ±2°F/1°C.
Resolution: 0.1°/1° resolution.
- Digital display with 0.1° resolution.
- Dimensions: 5.9 x 0.8 x 0.7" (150 x 20 x 18mm).
- Weight: 0.7oz (20g).
- Complete with one LR44 button battery and protective cover.



39272

39272, Pocket Fold-up Thermometer, -50°C to +300°C

4.5" stainless steel probe adjusts from 0 to 180° with detents at 45°, 90°, 135° and 180°.

- Wide °C/°F switchable temperature range -58 to 572°F (-50 to 300°C).
- Fold up probe 4.5" (114mm) adjusts from 0° to 180° for best measurement and viewing angle.
- Large LCD display with 0.1 degree resolution.
- Auto power off after 1 hour or when probe is placed in stored position.
- Data Hold and Min/Max recordings.
- Fast response time — one second under normal conditions.
- Complete with wrist strap and one AAA battery.

Full Range Of Thermometers: www.mrclab.com

Model	39272
Range	-58 to 572°F (-50 to 300°C)
Accuracy	±1.8°F (-22 to 482°F); ±1°C (-30 to 250°C)
Resolution	0.1°/1° resolution
Dimensions	6.1 x 2.0 x 0.8" (154 x 50 x 20mm)
Weight	2.5oz (71g)



392050

392050, Dual Stem Thermometer

Choice of Flat Surface or Penetration Stem Thermometer, For use in liquids, semi-solids and solids.

- Temperature range: -58 to 302°F (-50 to 150°C).
- Basic Accuracy of 2°F or 1°C.
- 0.1° resolution to 199.9; 1° over 200°.
- ON/OFF switch saves battery life.
- Auto power off after 1hr of non use.
- Dimensions: 5" (130mm) stainless steel probe with 2" (54mm) dial.
- Weight: 1oz (29g); 1.3oz (30g) with protective sheath.
- Includes 1.5V button battery, and plastic protective sheath with pocket clip.