



Laboratory Equipment Manufacturer
www.mrclab.com



A possible selection table of weights for calibration of weighing machines

	Resolution							
Capacity	100g	10g	1g	100mg	10mg	1mg	0.1mg	<0.1mg
Up to 50g		M3	M3	M3	M2	F2	E2	E1
Up to 100g	M3	M3	M3	M3	M1	F1	E1	E1
Up to 500g	M3	M3	M3	M2	F2	E2		
Up to 1kg	M3	M3	M3	M1	F1	E1		
Up to 5kg	M3	M3	M2	F2	E2			
Up to 10kg	M3	M3	M1	F1	E1			
Up to 50kg	M3	M2	F2	E2				
Up to 100kg	M3	M1	F1					
Up to 500kg	M2	F1	E2					

Note: This table should be interpreted in connection with 4.2.2 & 4.2.4 of the text