



TEPS-2



[www.mrclab.com](http://www.mrclab.com)

### TEPS-2, DIGITAL THERMOSTATIC HEATING SYSTEM MADE OF STAINLESS STEEL

- ▶ Length 206mm, width 112mm, height 236mm
- ▶ Weight 2.5kg
- ▶ Stainless steel case
- ▶ Input voltage 220V 50Hz, Integrated socket with 6amps fuse & circuit breaker
- ▶ PID controller for setting desired temperature
- ▶ 1000W heater power
- ▶ 35W 1440rpm engine with unique fluid mixing propeller
- ▶ Safety thermostat for extinguishing in case of fluid.