




K Series

# DSK700





K Series

**DSK700**

<b>Optical system</b>	Common Main Objective (CMO) Infinity optical system
<b>Observation tube</b>	Dual view Binocular head
<b>Inclination</b>	45° inclined
<b>Interpupillary distance</b>	54-76mm
<b>Diopter adjustment</b>	On both tubes, +/- 5 diopter
<b>Eyepieces</b>	Widefield WF10X/23mm
<b>Objectives system</b>	Zoom, ratio 5.2:1
<b>Objectives</b>	0.6X~3.1X
<b>Working distance</b>	89mm
<b>Stand type</b>	Pole type stand
<b>Head holder</b>	Optical bridge with light pointer
<b>Focus mechanism</b>	Coarse focusing system with tension adjustment
<b>Transformer</b>	External
<b>Power supply</b>	220-240V (CE)
<b>Accessories included</b>	Dust cover